Instructions for use

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**Syntesis<sup>®</sup>**

**phibo<sup>®</sup>**

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# Syntesis®

## PHIBO® SYNTESIS® CUSTOMISED ABUTMENTS

IMPORTANT INFORMATION. PLEASE READ CAREFULLY BEFORE USING PRODUCT

### INSTRUCTIONS FOR USE

#### Syntesis® customised abutments and custom anatomic abutments

Phibo®, Syntesis®, and Duplit®, are registered trademarks of Phibo Dental Solutions, S.L

This document contains information for the use of Phibo® prosthodontic components for dental implant systems and for prosthodontic components from other implant systems included in the Phibo® connection library.

#### PRODUCT DESCRIPTION

Customised abutments: customised abutment on implant made using CAD-CAM technology designed for use in cemented prosthetic rehabilitation. Materials: Titanium, Golden Titanium, 4-colour Zirconia, PMMA (Polymethyl methacrylate).

Custom anatomic abutments: ceramisable custom anatomic abutments (except PMMA) made using CAD-CAM technology, designed for use in screw-fitted prosthetic rehabilitation. Materials: Cobalt Chrome, 4-colour Zirconia, PMMA.

Syntesis® abutments allow for individualised customisation of function and aesthetics. Syntesis® abutments are designed by a Phibo® prosthetics specialist or by a laboratory that has been audited and validated by Phibo®. Design specifications are sent to Phibo®, where the customised abutment is manufactured.

**MODELS\***

ZR00XPP001CE	Customised abutment in ZR00 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
ZR01XPP001CE	Customised abutment in ZR01 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
ZR02XPP001CE	Customised abutment in ZR02 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
ZR03XPP001CE	Customised abutment in ZR03 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
ZR00XPD001CE	Custom Duplit in ZR00 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
ZR01XPD001CE	Custom Duplit in ZR01 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
ZR02XPD001CE	Custom Duplit in ZR02 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
ZR03XPD001CE	Custom Duplit in ZR03 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
TIXXXPP001CE	Customised abutment in Titanium connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
TIXXXPD001CE	Custom laboratory Duplit in Titanium connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
TIAUXPP001CE	Customised abutment in Golden Titanium connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
TIAUXPD001CE	Custom laboratory Duplit in Golden Titanium connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
PMA1XPP001CE	Customised abutment in A1 colour PMMA connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
PMB2XPP001CE	Customised abutment in B2 colour PMMA connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
PMA3XPP001CE	Customised abutment in A3 colour PMMA connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
PMA35PP001CE	Customised abutment in A3.5 colour PMMA connection PHIBO TSH / BNT (DI) SERIES 2, for cementing
ZR00XPP001AX	Customised abutment in ZR00 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, anatomic
ZR01XPP001AX	Customised abutment in ZR01 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, anatomic

ZR02XPP001AX	Customised abutment in ZR02 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, anatomic
ZR03XPP001AX	Customised abutment in ZR03 colour Zirconia connection PHIBO TSH / BNT (DI) SERIES 2, anatomic
PMA1XPP001AX	Customised abutment in A1 colour PMMA connection PHIBO TSH / BNT (DI) SERIES 2, anatomic
PMB2XPP001AX	Customised abutment in B2 colour PMMA connection PHIBO TSH / BNT (DI) SERIES 2, anatomic
PMA3XPP001AX	Customised abutment in A3 colour PMMA connection PHIBO TSH / BNT (DI) SERIES 2, anatomic
PMA35PP001AX	Customised abutment in A3.5 colour PMMA connection PHIBO TSH / BNT (DI) SERIES 2, anatomic
CCXXPP001AX	Customised abutment in Cobalt Chrome connection PHIBO TSH / BNT (DI) SERIES 2, anatomic

\*References vary depending on connection. An example of each abutment is shown with the Phibo® TSH/BNT series 2 direct to implant connection.

Thanks to the Phibo® connections library, we can also manufacture abutments for other implant systems.

**INDICATIONS**

In general, Syntesis® customised abutments are indicated to adapt or modify the transgingival contour and the emergency profile of the artificial restoration.

Syntesis® customised abutments in Zirconia, Titanium and Golden Titanium are indicated in unit, partial and total cemented restorations in all oral positions. These abutments are NOT INDICATED for direct ceramic coating.

Syntesis® customised abutments in PMMA are indicated in provisional cemented unit restorations (for no longer than 6 months). Specifically indicated in cases where prior modelling of soft tissue is required, before a definitive restoration device is fitted.

Syntesis® custom anatomic abutments in Zirconia and Cobalt Chrome are indicated in screw-fitted unit restorations for all oral positions where non-articulated continuity anatomic forms are required, depending on structural status and the health of the periodontal tissue, along with the patient's personal and clinical profile. These abutments are indicated for direct ceramic coating.

Syntesis® custom anatomic abutments in PMMA are indicated in provisional screw-fitted unit restorations (for no longer than 6 months), specifically when prior modelling of soft tissue is required, before fitting the definitive restoration.

Duplicate abutments ("Duplit" references) are indicated for use in dental laboratories for complete rehabilitation planning. They should not be used in the mouth and are not indicated for direct ceramic coating.

Syntesis® customised abutments are made for each patient. Accordingly, they are single-use products. Reuse of single-use products may involve their deterioration, at the risk of tissue infection, surgical or prosthodontic failure and/or deterioration of the patient's health.

After clinical and radiological diagnosis and the use of study models, implant treatment must be planned on the basis of the type of prosthodontic rehabilitation required by each patient, with a view to ensuring the success of treatment and expectations.

#### **CONTRAINDICATIONS**

In general, Syntesis® custom implants cannot be prescribed without prior prosthodontic and implant study and planning based on proper clinical diagnosis.

Custom anatomic abutments are contraindicated for screw-fitted unit crowns when the volume of the restoration and its height in relation to the antagonist exceed the patient's standard natural dimensions. In any case, a risk-benefit analysis is to be performed in relation to the implant to be used and the patient's anatomic conditions.

Custom anatomic abutments for screw-fitted unit crowns in Cobalt Chrome are contraindicated whenever the patient presents relative or absolute clinical contraindications requiring the use of other types of customised abutments:

- Relative contraindications: oral hygiene, periodontal pathologies, among others.
- Absolute contraindications: radiotherapy, treatment with immunosuppressant drugs, among others.

In the case of patients presenting local or systemic factors which could affect the restoration process, the prescribing clinician shall evaluate the potential benefits and risks of treatment.

#### **STORAGE AND HANDLING**

Phibo® Syntesis® products should be stored at between +10 and +40°C in a clean, dry place, protected from adverse conditions.

#### **RESPONSIBILITIES AND WARRANTY**

Syntesis® abutments are part of an overall concept and should only be used in conjunction with associated original products, in accordance with Phibo® instructions and guidelines and those of the manufacturers of other implant platforms compatible with Syntesis® (see the list of compatible implants at [www.phibo.com](http://www.phibo.com))

The use of third-party products in conjunction with Phibo® prosthetic components will invalidate all warranties and other obligations, whether explicit or implicit, provided by Phibo®. Likewise, the use of products manufactured by third parties other than those contained in the list of Phibo® Syntesis®-compatible implant systems in conjunction with prosthetic components from said manufacturers, will invalidate all warranties and other obligations, whether implicit or explicit, provided by Phibo®.

Phibo® Syntesis® customised abutments may be prescribed and fitted only by properly-trained professionals. In order for warranties to be valid, the clinician responsible for the treatment shall provide all information required by Phibo®.

Clinicians prescribing Phibo® prosthetic components will be responsible for establishing whether or not the product is indicated for each patient and each set of circumstances. Phibo® accepts no liability, either implicit or explicit, and will not be responsible for any damage, whether direct, indirect, punitive or of any other nature, produced by or arising in relation to errors of judgement or praxis by clinicians in the use of Phibo® products. Clinicians are also obliged to keep up to date with all the latest developments and applications regarding Phibo® products. If in doubt, please contact Phibo®. As the processing and surgical application of this product are under the control of the clinician, they are also the responsibility of the clinician. We accept no responsibility for damages arising from this cause.

#### PRECAUTIONS DURING THE PROCEDURE

Close collaboration between surgeon, clinician and laboratory technician are essential for the success of implant treatment. It is of particular importance to achieve the right balance of forces by means of the passive adaptation and adjustment of the bridge to the implant abutments, adjusting the occlusion of the opposing jaw and avoiding excessive transversal forces, particularly in cases of immediate loading.

#### CLINICAL PROCEDURE\*

##### Imprinting

1. Remove the provisional prosthesis or healing plug.
2. Insert the liner or open- or closed-tray imprint transfer on the implant or intermediate abutment.
3. Take the imprint, combining light silicon (high capacity for emergency profile copying) and heavy silicon. Remember to take the imprint of the antagonist, using putty or heavy silicon, and the bite.
4. If you have used a closed-tray imprint transfer, remove the transfers and place them correctly in the imprint, ensuring that the flat sides are in the right position to avoid rotation.
5. Whether you have used open-tray or closed-tray imprint transfers or liners, send the imprints to Phibo®, following the instructions given at [www.syntesisorder.com](http://www.syntesisorder.com).
6. Reinsert the provisional prosthesis or healing plug.

##### Fitting the abutment

7. To check the fit of the customised abutment, remove the provisional prosthesis or healing plug and insert the customised abutment, either on the implant or on an intermediate abutment.
8. After checking the fit, remove the abutment and send it to the lab for the definitive prosthesis to be made.

\*NB: There are two options here: the abutment may be connected directly to the implant or to an intermediate abutment.

Phibo® recommends that the definitive prosthesis also be manufactured using CAD-CAM technology, in order to guarantee the quality of the rehabilitation process.

#### DISINFECTION AND STERILISATION

Syntesis® prosthodontic components are not supplied sterile. Do not use products if packaging has been damaged or opened.

Syntesis® prosthodontic components must be cleaned, disinfected and sterilised prior to use. Syntesis® products should be disinfected using specific fast-acting products, with a pH of between 5 and 9, offering broad anti-viral and anti-fungal action. Follow manufacturer's instructions for use. Do not use peroxide, oxidising acids (nitric acid, sulphuric acid, oxalic acid, etc.) or products with a high chlorine content. After disinfection, rinse well in water and dry. Metal products (Titanium and Cobalt Chrome) should be sterilised in steam autoclaves, using a sterilising cycle at 134°C, for at least 6 minutes. For Zirconia abutments, the autoclave sterilisation cycle is as follows: at least 3 minutes and at most 6 minutes at 134°C. PMMA: due to its physical and chemical properties, the recommended sterilisation process is 15 minutes in autoclave at 121°C.

In general, we recommend disinfection with a high-level cold disinfectant. Follow the manufacturer's instructions for high-level disinfection. In the case of PMMA, we suggest disinfection for 10 minutes with

2% sodium hypochlorite. Do not remove the sterilised product until after completion of the drying cycle. We recommend using a sterilisation control, recording the date and expiry date, in addition to performing periodical controls of the sterilisation process using biological indicators.

**IMPORTANT:** Do not sterilise plastic products in the autoclave. Do not use dry heat sterilisers, as they may damage metal and plastic products. Plastic products may be sterilised before use by cold disinfection using chemical disinfectants.

#### **PRECAUTIONS**

##### **PROSTHODONTIC REHABILITATION**

Prosthetic procedures describe in detail the precautions to be observed during treatment. The rehabilitation process and the prosthesis should be designed before implants are inserted.

#### **WARNING**

Syntesis® abutments are custom-designed for maximum functionality and aesthetics and do not require modification. Subsequent modifications may affect their mechanical properties and the properties of the materials used. Modifications to abutments may invalidate product warranties.

Use only the Syntesis® screws provided with the abutment. This screw may not be replaced by a screw not manufactured by Phibo®.

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